

H.pylori Ab TEST



Infection with the gram-negative *Helicobacter pylori* bacteria elicits an immune response in the human body including the production of specific anti- *H.pylori* antibodies in blood and other body fluids.

The antibodies aid the body in the fight against the infection and can be used for the detection of past and possibly active *H.pylori* infections.

INTENDED USE

Turklab's *H.pylori* Ab Test is a rapid, immunochromatographic assay for the qualitative detection of the antibodies against the *H.pylori* antigens in human whole blood, serum or plasma.

- If antibodies are present, it may mean that they are there to fight *H. pylori* bacteria.

PERFORMANCE EVALUATION

- Sensitivity :98,7%
- Specificity :93,2 %



PRODUCT	TEST TYPE	SAMPLE TYPE	READING TIME	KIT SIZE	KIT COMPONENTS
H.Pylori Ab	CASSETTE	WB/S/P	10 Min.	40	Test cassette + Pipette + Dilution Buffer + Instructions For Use