

## Rotavirus Adenovirus TEST



Rotavirus and Adenovirus are two infectious agents that cause mild to severe diseases.

These viral pathogens have been isolated throughout the world where the infections are most frequently seen in children less than two years of age but have been found in patients of all ages. Rapid and accurate diagnosis of gastroenteritis due to Adenovirus and rotavirus helps establish the etiology of gastroenteritis and related patient management.

### INTENDED USE

Turklab's *Rotavirüs Adenovirus Test* is a rapid, immunochromatographic assay for the qualitative detection of *Rotavirus* and *Adenovirus* antigens in human feces samples.

- Clean and reliable use.
- Ease of use with the sample collection tube with prefilled diluent.

### PERFORMANCE EVALUATION

#### Rotavirus:

- Sensitivity :99,9%
- Specificity :98,4%

#### Adenovirus:

- Sensitivity :99,9%
- Specificity :99,0%



PRODUCT	TEST TYPE	SAMPLE TYPE	READING TIME	KIT SIZE	KIT COMPONENTS
ROTAVIRUS-ADENOVIRUS	CASSETTE	FECES	10 Min.	40	Test cassette + pipette + sample collection tube with extraction buffer + instructions for Use